visocor[®] OM50

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(DE-2 - DE-35)

Instructions for Use

(EN-36 - EN-69)

www.visocor.de • *Besser Messen.*

UEBE[®]
Germany

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Table of contents

A	Mode of operation	38
B	Safety instructions	
1.	Important instructions for patients	39
2.	Important technical details	40
C	Operating the unit	
1.	Unit description	42
2.	Displays	43
3.	Important instructions for use	44
4.	Initial operation of the unit	45
5.	Inserting/replacing batteries	45
6.	Activating and setting the date/time	46
7.	Attaching the cuff	48
8.	Measuring blood pressure	49
9.	Irregular pulse waves	50
10.	Using the memory	51

Table of contents

D	What you should know about blood pressure	
1.	Systolic and diastolic blood pressure values	53
2.	Why you measure different values	53
3.	Why you should measure blood pressure regularly	53
4.	What are normal values for blood pressure?	54
E	Technical information	
1.	Failure and error messages	55
2.	Customer service	57
3.	Technical data	58
4.	Original spare parts and accessories	59
5.	Applicable standards	60
6.	Technical inspection / calibration check	61
F	General provisions	
1.	Explanation of symbols	62
2.	Disposal	63
3.	Electromagnetic compatibility (EMC)	63
G	Maintaining the unit	68
H	Warranty	69

A Mode of operation

Thank you for choosing the visocor® OM50 upper arm blood pressure monitor (also referred to in the following as the unit). This unit is recommended for patients with unstable blood pressure, for measuring their blood pressure at home and for therapy support.

visocor® OM50 uses the oscillometric method for measuring blood pressure and pulse rate.

Using the pressure variations caused by the pulse at falling cuff pressure, the micro-computer calculates the values systole, diastole, and pulse. Afterwards, the cuff is deflated completely.

A measurement value memory stores the most recent 60 measurement results in each case.

These instructions are intended to help you use the unit reliably and effectively. Keep the instructions together with the unit at all times. The unit must be used in accordance with the methods detailed in these instructions and may not be used for other purposes. Read these instructions carefully before using the unit.

Safety instructions B

1. Important instructions for patients

- The unit is designed for non-invasive measurement of the systolic and diastolic blood pressure on the upper arm, as well as measurement of the pulse rate of adults, i.e. 15 years of age and older. Blood pressure measurements on children require specialist knowledge! Please consult your doctor if you wish to measure the blood pressure of a child. Do not under any circumstances use the unit on a baby/infant.
- The cuff has been specially developed for this unit and must not be used with other units.
- The measurement results of automatically measuring blood pressure monitors may be falsified by pregnancy, irregular heartbeat or arteriosclerosis. Measure your own blood pressure in consultation with your doctor.
- Do not under any circumstances place the cuff on or over any critical point, e.g. wound, aneurysm, etc. Risk of injuries! The supply via an intravascular access (infusion) might be interrupted.
- The unit contains small parts which could be swallowed by children. It should therefore be kept out of the reach of children at all times.
- If you have had a mastectomy, do not measure blood pressure on the side of the body where the breast has been removed together with the lymph nodes from the armpit.
- Measuring your own blood pressure does not constitute treatment. On no account modify of your own accord the dosage of drugs/medication prescribed by your doctor.

B Safety instructions

- Please refer to the chapter "Important instructions for use" (Page EN-44) before conducting any measurements yourself.

2. Important technical details


- A consistently good power supply to the unit is essential for trouble-free blood pressure measurements.
 - Please use long-life alkaline batteries only (LR6).
 - When replacing batteries, always replace all the batteries at the same time.
 - You will need 4 x 1.5 V AA/LR6 batteries. Rechargeable batteries have a voltage of only 1.2 V and are thus not suitable.
 - When using the unit with a mains adapter, please use only the type A1 mains adapter specially tested for medical devices.
 - If the unit is not used for an extended period of time, please remove the batteries. Essentially, all batteries can leak.
- The unit may only be operated with original spare parts. The warranty will be invalidated if the unit is damaged by non-approved accessories!
- Using this unit near mobile phones, microwave ovens or other devices with strong electromagnetic fields may lead to malfunctions. See also EMC description from Page EN-63.
- The display of the pulse frequency is not suitable for checking the frequency of cardiac pacemakers. Cardiac pacemakers and blood pressure monitors do not influence each other in their mode of operation.

Safety instructions B

- Never open or modify the unit or the cuff - this is a medical device and may only be opened by authorised specialists. (Exception: replacing batteries). If the unit has been opened, it must be subjected to a metrological inspection by an authorised institution.
- The cuff may be inflated on the arm only.
- Please comply with the specified operating conditions for measurement. See Technical data, Page EN-58.
- The inflating and measuring procedure can be stopped by pressing the Start/Stop button. In this case, the unit stops the inflation procedure and deflates the cuff.



C Operating the unit

1. Unit description

C Operating the unit

2. Displays

Failure and error messages (Page EN-55)

C Operating the unit

3. Important instructions for use

- Refrain from drinking alcoholic or caffeinated beverages and smoking at least one hour before measuring.
- Please rest at least 5 minutes before measuring. Depending on the previous degree of stress/exertion, this can even require up to one hour.
- Expose the upper arm. Clothing must not impair the blood flow in or from the arm, as this influences the blood pressure at the measurement point.
- Body posture must be relaxed:
 - For this purpose, sit at a table (if possible the height of a dining table, not a coffee table!).
 - Rest your back against the backrest of the chair.
 - Rest your forearm completely flat, with the palm of your hand facing upwards.
 - Place your feet on the floor and do not cross your legs.
- There should be no irregular heartbeats during the measurement! Irregular movements, vibrations, speaking, and breathing heavily will also affect the measurement.
- It is absolutely essential to sit still and remain quiet. Irregular movements, vibrations (shaking), talking and breathing heavily will affect and may even falsify the measurement. Look out for irregular pulse signals in the display; if necessary, repeat the measurement under more favourable conditions.
- Unusual measurement results are possible with all automatically measuring blood pressure measurement systems from time to time. Check yourself: did you comply with the instructions for use above? If necessary, repeat the measurement after allowing the blood circulation in the arm to recuperate briefly; rest for about 3-5 minutes for this purpose.
Tip: Remain seated, your monitor will switch off automatically about 3 minutes after a measurement. Afterwards, we recommend that you repeat the measurement.

Operating the unit C

- Blood pressure is not a fixed value. It may fluctuate upwards or downwards by more than 20 mmHg on patients within a few minutes.

4. Initial operation of the unit

Insert the supplied batteries into the unit, if these have not been inserted already.


If the unit is to be operated from the mains power supply, the cable plug of the mains adapter (special accessory) must be inserted into the connection socket on the back of the unit. Please use the UEBE type A1 mains adapter only. The batteries are switched off automatically.

5. Inserting/replacing batteries

- Opening the battery compartment
Remove the cover of the battery compartment on the underside of the unit.
- Inserting batteries
Remove the old batteries from the unit and insert the new ones. Please observe proper polarity (marking in battery compartment).
- Closing the battery compartment
Close the battery compartment by clipping the battery cover back into the unit.

C Operating the unit

When the display shows the „Flat battery“ symbol for the first time, the unit will still have enough charge for about 30 measurements. Please change the batteries during this period.

6. Activating and setting the date/time


The unit has a date/time function which you can activate where necessary. This function is deactivated when the unit is delivered.

Activating the date/time function

To activate the date/time function, press the Start/Stop button and then immediately, while the full display is visible, the Memory button. This function can be deactivated at any time in the same way.

Setting the date/time

The four-digit number representing the year appears in the display first. Press the Memory button to set the year. Press the Start/Stop button to confirm the selected year and proceed to setting the date.

C Operating the unit

Set the month with the Memory button and confirm this with the Start/Stop button. Set the day, hours and minutes in the same way.



Changing the date/time

To change the date/time function, take out a battery and wait until the display goes blank. Now you can reactivate the date/time function and reset the values.

Deactivating the date/time function

To deactivate the date/time function, press the Start/Stop button and then immediately the Memory button.

When the date/time function is activated, the time is also shown in the display when the unit is switched off.

C Operating the unit

7. Attaching the cuff

- Expose the upper arm.
- Push the cuff onto the arm (Fig. 1) until the lower edge of the cuff is 2-3 cm above the crook of the arm (Fig. 2).
- When using the unit on the left arm, the air hose runs to the unit through the middle of the crook of the arm so that the 4 cm long artery marking is automatically located centrally above the pulse position (Fig. 3).
- When using the unit on the right arm, the cuff must be rotated to the left until the artery marking is on the pulse position. The air hose then runs along the inside of the upper arm.
- The cuff should be tightened to such an extent that 2 fingers can still be inserted between arm and cuff. Tip:
 - Bend your arm slightly.
 - The upper arm muscle is tensioned slightly.
 - This slightly increases the circumference of the arm.
- Now pull the free cuff end tight and close the Velcro fastener (Fig. 4)


Fig. 1


Fig. 2


Fig. 3



Fig. 4

C Operating the unit

- Check if the measurement arrow on the cuff is inside the „marking for arm circumference“ at the edge of the cuff.
- Lay your extended and relaxed arm with the cuff on the dining table and do not move it during the measurement; do not talk. The palm should be facing upwards.
- Connect the cuff connector to the cuff socket on the left-hand side of the unit. Make sure the connector is fully inserted into the unit.


8. Measuring blood pressure

- Switch on the unit by pressing the Start/Stop button. The following display appears:
- After the unit has completed its calibration against the ambient air pressure, the cuff is automatically inflated to approx. 180 mmHg. If that cuff pressure is not sufficient for measuring the blood pressure level, the device automatically increases the cuff pressure in steps of 30 mmHg until a sufficient cuff pressure has been reached.
- For an inflation pressure exceeding 180 mmHg the inflation process can also be controlled manually: Keep the Start/Stop button pressed immediately after the inflation has started and release it when the cuff pressure has reached 40 mmHg above the expected maximum systolic blood pressure. The inflation stops with releasing the button.


C Operating the unit

- The actual measurement takes place with the start of deflation. The "♥" symbol in the display flashes until the pulse frequency is displayed.


- The measurement being completed is indicated by a long beep. The determined values for systole, diastole, and pulse are shown on the display.
- The unit switches off automatically after about 3 minutes.

9. Irregular pulse waves

If the symbol for irregular pulse waves flashes after a measurement, the unit has recorded irregular pulses during the measurement. This can be set off by an irregular heartbeat (arrhythmia), disturbance caused by movement, talking or even breathing in deeply. The symbol is stored with the respective measurement.


Operating the unit C

10. Using the memory

Measured results are automatically stored in the memory. The memory can store up to 60 results and the average value. When more than 60 measured values have been stored, the oldest value (No. 60) is deleted to allow the latest value (No. 1) to be recorded.


Recalling data

Press the Memory button to call up data. The average value of the stored results is displayed with „A“.

When the Memory button is pressed again, the value measured last (No. 1) appears.

Press the Memory button repeatedly to call up further measured values.

When the date/time function is activated, the measured values are displayed alternately with memory number, date and time.

The unit switches itself off automatically 30 seconds after the memory is called up.

C Operating the unit

Deleting data

To delete individual values, press the Memory button several times to obtain the measured value required.

Then press the Memory button again and hold it down for 8-10 seconds until the value first flashes and then disappears.

To delete the entire memory, recall the relevant average value display (A) and press the MEM button again, holding it down until the mean value first flashes and then disappears.

What you should know about blood pressure **D**

1. Systolic and diastolic blood pressure values

The cardiovascular system has the important function of supplying all organs and tissues in the body with sufficient amounts of blood and of transporting metabolites. To do so, the heart contracts and expands at a regular rate of about 60 to 80 times per minute. The pressure of the flowing blood on the artery walls caused by the heart contracting is termed systolic. The pressure in the ensuing relaxation phase, when the heart refills with blood, is termed diastolic. During daily measurement you determine both values.

2. Why you measure different values

Our blood pressure responds to internal and external influences like a sensitive measuring instrument. It can be affected by even slight changes. This explains why values measured by a doctor or pharmacist are higher than those measured at home in the environment you are used to. Changes in the weather, climate changes, or physical or psychological stress can have effects as well.

3. Why you should measure blood pressure regularly

Even the time of day has an influence on your blood pressure. During the day the values are generally higher than during the periods of rest at night. One-off and irregular measurements therefore say little about your actual blood pressure. A reliable assessment is possible only when measurements are taken regularly. Discuss the measurement values with your doctor.

D What you should know about blood pressure

4. What are normal values for blood pressure?

The World Health Organisation (WHO) established the following limits in mmHg (millimetre hydrargyrum) to classify the blood pressure values.

WHO 2003	Systolic pressure = Upper value mmHg**	Diastolic pressure = Lower value mmHg**
High pressure	from 140*	from 90*
Limit	120 to 139	80 to 89
Normal value	below 120	below 80

* It is sufficient for one of the values to be raised.

** millimetres of mercury

Technical information E

1. Failure and error messages

Failure encountered	Possible cause	Corrective action
Display Err - 300	Excessive cuff pressure. Cuff inflated to maximum due to movement of the arm or body.	- Repeat measurement - Do not move arm - Do not talk
Display Err - 1	Measurement error. Measurement could not be carried out due to movement of the arm or body.	- Repeat measurement - Do not move arm - Do not talk
Display Err - 2	Inflation error, pressure built up too quickly or too slowly.	Check that air connector is correctly fitted, repeat measurement.
	Fault was recognised.	- Repeat measurement - Do not move arm - Do not talk
Display Err - 3	Deflation rate too high or too low.	Call the Service Hotline if this occurs repeatedly.
	Irregular heartbeat, disruptive movements, shaking, wobbling (objects), breathing in deeply, etc.	Repeat measurement after 3-5 minutes rest. Possible factors (depending on the severity of the arrhythmia) affecting the measurement results must be discussed with your doctor.

E Technical information

Failure encountered	Possible cause	Corrective action
Measured values too high	Too little rest before the measurement.	Repeat measurement after a break of approx. 3-5 mins. - Do not move arm - Do not talk
Unusual measured values	Movement or talking during measurement, resting time not observed, feet possibly crossed, smoking or drinking coffee.	Check conditions and repeat measurement. Observe instructions for use P. EN-44.
Display 0 0	Start/Stop button accidentally pressed while battery was being changed.	Switch unit off and then on again with Start/Stop button; if necessary.
Measurement was interrupted and cuff first deflated and then inflated	Unit detects a fault or that diastolic value is too low.	Unit shows no fault. Repeat measurement.
	Movement during measurement	If necessary, interrupt measurement, measure again after a break of 5 minutes.
No display after unit is switched on	Batteries inserted incorrectly.	Check position of batteries.
	Batteries flat.	Change batteries.
	Battery contacts dirty.	Clean battery contacts.
No display after Memory button is pressed	No measured values stored.	Take a new measurement.
Measurement interrupted	Batteries flat	Change batteries.

Technical information E

Failure encountered	Possible cause	Corrective action
Measured values extremely high or low	Incorrect measurement position.	Repeat measurement, ensuring correct body posture. Do not move your arm and do not talk.
Date/time not shown in display	- Function accidentally deactivated - Function was not reactivated after battery replacement	Activate date and time as described on P. EN-46.

2. Customer service

The unit may only be repaired by the manufacturer or by an expressly authorised organisation. Please contact:

UEBE Medical GmbH
 Zum Ottersberg 9
 97877 Wertheim, Germany
 Tel.: +49 (0) 9342/924040
 Fax: +49 (0) 9342/924080
 E-mail: info@uebe.com
 Internet: www.uebe.com

E Technical information

3. Technical data

Unit type:	Digital automatic unit with electric pump for measuring blood pressure on upper arm
Dimensions:	L = 100 mm x W = 152 mm x H = 61 mm
Weight:	295 g not including batteries
Display:	LCD display (liquid crystal display) for measured values and check displays
Memory:	60 measured values (stored automatically) and average value (A)
Measurement procedure:	Oscillometric determination of systole, diastole and pulse
Reference procedure of clinical testing:	Auscultatory measurement
Inflating pressure:	Initially up to approx. 180 mmHg, then in stages of 30 mmHg
Pressure display range:	0-300 mmHg
Measurement range	Systolic: 50-250 mmHg Diastolic: 40-150 mmHg Pulse: 40-160 pulses/min
Error limits:	Blood pressure measurement: corresponds to EN 1060 Part 3 Pressure measurement: ± 3 mmHg Pulse measurement: $\pm 5\%$

Technical information E

Serial number:	The unit bears a serial number SN which provides clear identification.
Power supply:	Battery type: 4 x 1.5 V AA alkali manganese (LR 6) or lithium (FR 6) round cells, service life: more than 800 measurements in 2 years
	Optionally stabilised type A1 mains adapter, output 6 VDC, minimum 600 mA
Cuff:	Cuff VM5 (for arm circumference 22-32cm) or optionally cuff VL5 (for arm circumference 32-42 cm)
Operating conditions:	Ambient temperature 10 to 40 °C, relative humidity up to 85 %, non-condensing
Storage and transport conditions:	Ambient temperature -20 to +50 °C, relative humidity up to 85 %, non-condensing
Air deflation valve:	Electronically controlled
Automatic switch-off	approx. 3 minutes after end of measurement
IP code:	IPX0 not protected

4. Original spare parts and accessories:

- Cuff 22-32 cm type VM5
Part. no. 2505001
PZN 9101599
- Cuff 32-42 cm type VL5
Part. no. 2505002
PZN 9101607
- Mains adapter type A1
Part. no. 2501020
PZN 0179482

Subject to technical modifications.

E Technical information

5. Applicable standards

- EN 1060-1 : 1995 + A2 : 2009 Non-invasive sphygmomanometers - Part 1: General requirements
- EN 1060-3 : 1997 + A2 : 2009 Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems
- IEC 60601-1:2005 Medical electrical equipment - Part 1: General requirements for basic safety
- IEC 60601-1-2 : 2007 Medical electrical equipment
Electromagnetic compatibility
- ISO 81060-2:2009 Non-invasive sphygmomanometers - Part 2: Clinical validation of automated version
- DIN EN ISO 10993-1:2009 Biological evaluation of medical devices - Part 1: Evaluation and testing

Manufacturer: UEBE Medical GmbH
Zum Ottersberg 9
97877 Wertheim
Germany

Technical information **E**

6. Technical inspection / calibration check

Basically, a metrological inspection is recommended at intervals of 2 years. However, professional users in Germany are obligated to comply with the aforementioned according to „Regulation for Operators of Medical Devices“.

This can be implemented either by UEBE Medical GmbH, an authority responsible for metrology, or authorised maintenance services. For this, please observe your national provisions.

Upon request, responsible authorities or authorised maintenance services receive a „Test instruction for metrological inspection“ from the manufacturer.

Important: No modifications, e.g. opening the unit (except to replace the batteries) may be made to this unit without the manufacturer's prior permission.

F General provisions

1. Explanation of symbols

CE 0123 This product complies with the Council Directive 93/42/EC from 5 September 2007 regarding medical devices, which became effective on 21 March 2010 and bears the mark CE 0123 (TÜV SÜD Product Service GmbH). Units with CE marks are subject to quality inspections in accordance with this Directive and provide a higher level of accuracy than previous calibration.

Degree of protection against electric shock: TYPE BF

Please observe instructions for use

Storage and transportation conditions
Ambient temperature -20 to +50 °C

Relative humidity up to 85 %, non-condensing

Manufacturer

If using the optionally available mains adapter: Equipment class 2 (double-insulated)

General provisions F

2. Disposal

Batteries and technical appliances must not be disposed of with domestic waste, but should be handed in at the appropriate collection and disposal points.

3. Electromagnetic compatibility (EMC)

Technical description

The unit satisfies the EMC requirements of the international standard IEC60601-1-2. The requirements are satisfied under the conditions described in the tables below. The unit is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use. Portable and mobile HF communications equipment can affect the unit. Use of the unit in conjunction with non-approved accessories can affect the unit negatively and alter the electromagnetic compatibility. The unit should not be used directly adjacent to or between other electrical equipment.

F General provisions

Guidance and manufacturer's declaration - electromagnetic emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Complies	

General provisions F

Guidance and manufacturer's declaration - electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment. Note: Power plug type is 2-pin.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	(>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec.	(>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device is powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment..
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

F General provisions

Guidance and manufacturer's declaration - electromagnetic immunity			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF radio equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = [3.5/V1]\sqrt{P}$ $d = [3.5/E1]\sqrt{P}$ 80 MHz to 800 MHz $d = [7/E1]\sqrt{P}$ 800 MHz to 2.5 GHz
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, c should be less than the compliance level in each frequency range. d Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a: Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio stations and TV stations cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.			
b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

General provisions F

Recommended separation distances between portable and mobile RF communications equipment and the device.			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = [3.5/V1]\sqrt{P}$	80 MHz to 800 MHz $d = [3.5/E1]\sqrt{P}$	800 MHz to 2.5 GHz $d = [7/E1]\sqrt{P}$
0.01		0.12	0.23
0.1		0.37	0.73
1		1.2	2.3
10		3.7	7.3
100		11.7	23.3
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE1: At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

G Maintaining the unit

- The unit contains sensitive parts and must be protected against strong variations in temperature, air humidity, dust and direct sunlight.
- The unit is neither impact-resistant nor shock-proof. We recommend that you have the intactness and accuracy of the display checked after heavier falls or impacts.
- The unit is not waterproof.
- If the unit becomes dirty during use, use a soft, dry cloth to clean the unit. Do not use benzene, thinners or other strong solvents.
- Do not scrub or machine-wash the cuff. If the cuff gets dirty during use, use a synthetic cleaner and softly rub the surface.
- If the unit has been opened, it must be subjected to a metrological inspection conducted by an authorised organisation.
- Make sure that liquid cannot get into the air hose. Dry thoroughly.

Warranty **H**

The blood pressure monitor has been manufactured and tested with great care. However, in the unlikely event of a defect being detected after delivery, we provide warranty in accordance with the following terms and conditions:

During the warranty period of 2 years from the date of purchase we reserve the right either to repair any such defect at our expense upon return of the unit to our factory or to supply a perfect replacement unit.

Excluded from the warranty are parts subject to normal wear and tear as well as damage caused by non-compliance with the instructions for use, improper handling (e.g. unsuitable power sources, breakages, leaking batteries) and/or disassembly of the unit by the purchaser. Furthermore, no claims for damages against us are substantiated by the warranty.

Warranty claims can only be advanced in the warranty period and by presenting proof of purchase. In the event of a warranty claim, the unit must be sent to the following address together with the proof of purchase and a description of the complaint:

UEBE Medical GmbH
Service-Center
Zum Ottersberg 9
97877 Wertheim
Germany

In the case of defectiveness of the goods, the contractual rights of the purchaser to claim against the seller in accordance with § 437 German Civil Code are not limited by the warranty.

Please note:

In the event of a warranty claim it is essential to attach the proof of purchase.

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© Copyright 2011 UEBE Medical GmbH