BM 85

(D)	Blutdruckmessgerät Gebrauchsanweisung	2-18
(GB)	Blood pressure monitor Instructions for use	19-35
(F)	Tensiomètre Mode d'emploi	35-52
(E)	Tensiómetro Instrucciones de uso	52-69
(I)	Misuratore di pressione Istruzioni per l'uso	69-86
(TR)	Tansiyon ölçme cihazı Kullanım Kılavuzu	86-102
(RUS)	Прибор для измерения кровяного давления в плечевой артерии Инструкция по применению	102-122
(PL)	Ciśnieniomierz Instrukcja obsługi	121-139
	Electromagnetic Compatibility Information	138-141

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Dear Customer,

Thank you for choosing one of our products. Our name stands for high-quality, thoroughly tested products for applications in the areas of heat, weight, blood pressure, body temperature, pulse, gentle therapy, massage and air. Please read these instructions for use carefully and keep them for later use, be sure to make them accessible to other users and observe the notes they contain.

With kind regards,
Your Beurer team

1. Getting to know your instrument

Check that the packaging of the Beurer BM 85 blood pressure monitor has not been tampered with and make sure that all the required contents are present. Before use, ensure that there is no visible damage to the device or accessories and that all packaging material has been removed. If you have any doubts, do not use the device and contact your retailer or the specified Customer Services address.

The upper arm blood pressure monitor is used to carry out non-invasive measurement and monitoring of the arterial blood pressure values in adults.

This allows you to quickly and easily measure your blood pressure, save the measurements and display the development and average values of the measurements taken.

You are also warned of possible existing cardiac arrhythmia.

The recorded values are classified according to the WHO guidelines and evaluated graphically.

This blood pressure monitor also has a haemodynamic stability display, which is referred to as a resting indicator throughout these instructions for use. This shows whether you, and consequently your circulatory system, are sufficiently at rest when the blood pressure measurement is being taken and is therefore a more precise indicator of your resting blood pressure. Read more about this in section 6.

Store these instructions for use for future reference and make them accessible to other users.



2. Important notes

Signs and symbols

The following symbols are used in these instructions for use, on the packaging and on the type plate for the device and the accessories:

	Attention
	Note Note on important information
	Observe the instructions for use

	Application part, type BF
	Direct current
	Disposal in accordance with EC Directive WEEE (Waste Electrical and Electronic Equipment)
	Manufacturer
Storage 	Permissible storage temperature and humidity
Operating 	Permissible operating temperature and humidity
	Protect from moisture
SN	Serial number
	The CE labelling certifies that the product complies with the essential requirements of Directive 93/42/EEC on medical products.

Notes on use

- In order to ensure comparable values, always measure your blood pressure at the same time of day.
- Before every measurement, relax for about five minutes.

- If you want to perform several measurements on the same person, wait five minutes between each measurement.
- Do not take a measurement within 30 minutes of eating, drinking, smoking or exercising.
- Repeat the measurement if you are unsure of the measured value.
- The measurements taken by you are for your information only – they are no substitute for a medical examination! Discuss the measurements with your doctor and never base any medical decisions on them (e.g. medicines and their administration)!
- Do not use the blood pressure monitor on newborns or patients with preeclampsia. We recommend consulting a doctor before using the blood pressure monitor during pregnancy.
- Cardiovascular diseases may lead to incorrect measurements or have a detrimental effect on measurement accuracy. The same also applies to very low blood pressure, diabetes, circulatory disorders and arrhythmias as well as chills or shaking.
- The blood pressure monitor must not be used in connection with a high-frequency surgical unit.
- Only use the unit on people that have the specified upper arm measurement for the unit.
- Please note that when inflating, the functions of the limb in question may be impaired.
- During the blood pressure measurement, the blood circulation must not be stopped for an unnecessarily long time. If the device malfunctions remove the cuff from the arm.

- Avoid any mechanical restriction, compression or bending of the cuff line.
- Do not allow sustained pressure in the cuff or frequent measurements. The resulting restriction of the blood flow may cause injury.
- Make sure that the cuff is not placed on an arm in which the arteries or veins are undergoing medical treatment, e.g. intravascular access or intravascular therapy, or an arteriovenous (AV) shunt.
- Do not use the cuff on people who have undergone a mastectomy.
- Do not place the cuff over wounds as this may cause further injury.
- Please note that data transfer and data storage is only possible when your blood pressure monitor is supplied with power. As soon as the battery is empty, the blood pressure monitor loses the date and time setting.
- To conserve the battery, the blood pressure monitor switches off automatically if you do not press any buttons for 3 minutes.
- The device is only intended for the purpose described in these instructions for use. The manufacturer is not liable for damage resulting from improper or careless use.

Instructions for storage and maintenance

- The blood pressure monitor is made from precision and electronic components. The accuracy of the measurements and service life of the device depend on its careful handling:
 - Protect the device from impacts, humidity, dirt, marked temperature fluctuations and direct sunlight.
 - Do not drop the device.

- Do not use the device in the vicinity of strong electromagnetic fields and keep it away from radio systems or mobile telephones.
- Only use the cuff included with the delivery or original replacement parts. Otherwise incorrect measurements will be recorded.
- Do not press the buttons before the cuff is placed on the arm.

Notes on batteries

- Swallowing batteries can be extremely dangerous. Keep the batteries and products out of the reach of small children. If a battery has been swallowed, seek medical attention immediately.
- Batteries must not be taken apart, thrown into an open fire or short circuited.

The codes below are printed on batteries containing harmful substances:
 Pb: battery contains lead,
 Cd: battery contains cadmium,
 Hg: battery contains mercury.


Instructions for repairs and disposal

- Do not repair or adjust the device yourself. Proper operation can no longer be guaranteed in this case.
- Repairs must only be carried out by Customer Services or authorised suppliers.
- Do not open the device, except during disposal. In this case, remove the installed battery. Failure to comply will invalidate the warranty.

- Batteries must not be disposed of with household waste. Please dispose of batteries at the collection points intended for this purpose.
- Remove the battery before disposing of the device. Undo the four, round rubber covers on the rear of the device to remove the battery. Unscrew the housing. Remove the battery and dispose of it correctly.
- Please dispose of the device in accordance with EC Directive – WEEE (Waste Electrical and Electronic Equipment).
If you have any queries, please contact the appropriate local authorities.


3. Unit description

- Cuff
- Cuff line
- Cuff connector
- Cuff holder
- USB interface
- Display
- Memory buttons **M1/M2**
- START/STOP button **ⓘ**
- Resting indicator display
- WHO scale
- Connection for cuff connector (left-hand side)

Information on the display:

1. Time/date
2. Systolic pressure
3. Diastolic pressure
4. Calculated pulse value
5. Pulse symbol
6. Release air (arrow)
7. Number of memory space/memory display average value (M), morning (AM), evening (PM)
8. WHO classification
9. Battery level symbol
10. Cardiac arrhythmia symbol
11. User memory
12. Symbol for *Bluetooth*® transfer

USB interface

The blood pressure monitor also allows you to transfer your measured values to the PC.

To do this you require a commercially available USB cable (included in delivery) and the Beurer “HealthManager” PC software.

The software can be downloaded free of charge from the download area under Service at www.beurer.com.

System requirements for the Beurer “HealthManager” PC software.

1. Supported operating systems:
 - Windows XP SP3
 - Windows Vista SP1 or later
 - Windows 7
 - Windows 7 SP1
 - Windows 8
2. Supported architectures:
 - x86 (32 bit)
 - x64 (64 bit)
3. Hardware requirements:
 - Recommended: at least Pentium 1 GHz or faster with at least 1 GB RAM
 - Free memory on the primary partition of at least:
 - x86 – 600 MB
 - x64 – 1.5 GB
 - Graphic resolution from: 1024 x 768 pixels
 - USB port 1.0 or later

4. Preparing for the measurement

Completely charge up the blood pressure monitor before taking your first measurement:


- Using a PC:
Insert the USB cable into the USB interface on the device and directly connect it to your PC.
- Using the mains part
Insert the USB cable into the blood pressure monitor and connect it to mains supply using the enclosed USB mains part.

The  symbol is shown once the device is fully charged.

Charge the battery if the battery level symbol flashes. If the device's battery has been completely drained and you are unable to switch it on again, you must reset the date, time and *Bluetooth*[®] settings. In this process, any saved measurements are retained.

Set the hour format, date, time and *Bluetooth*[®] settings

The following section describes the functions and settings available on the blood pressure monitor.

It is essential to set the correct date and time. Otherwise, you will not be able to save your measured values correctly with a date and time and access them again later.

- If you press and hold the **M1** or **M2** memory button, you can set the values more quickly.

Press and hold the **START/STOP** button  for 5 seconds.

Hour format

The hour format now flashes on the display.

- Select the desired hour format using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

Date

The year flashes on the display.

- Select the year using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

The month flashes on the display.

- Select the month using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

The day flashes on the display.

- Select the current day using the **M1/M2** memory buttons and confirm with the **START/STOP** button .

- If you have set the 12h hour format, the month is displayed before the day.

The hours flash on the display.

- Select the current hours using the **M1/M2** memory buttons and confirm with the **START/STOP** button **(i)**.

The minutes flash on the display.

- Choose the current minutes using the **M1/M2** memory buttons and confirm with the **START/STOP** button **(i)**.

The *Bluetooth*® symbol flashes on the display.

- Use the **M1/M2** memory buttons to choose whether automatic *Bluetooth*® data transfer is activated (*Bluetooth*® symbol flashes) or deactivated (*Bluetooth*® symbol is not shown) and confirm with the **START/STOP** button **(i)**.

- (i)** *Bluetooth*® transfers will reduce the battery capacity.


5. Measuring blood pressure

Please ensure the unit is at room temperature before measuring.


The measurement can be performed on the left or right arm.

Attaching the cuff

Place the cuff on to the bare left upper arm. The circulation of the arm must not be hindered by tight clothing or similar.


The cuff must be placed on the upper arm so that the bottom edge is positioned 2–3 cm above the elbow and over the artery. The line points to the centre of the palm.

Now tighten the free end of the cuff, but make sure that it is not too tight around the arm and close the hook-and-loop fastener. The cuff should be fastened so that two fingers fit under the cuff.

Now insert the cuff line into the connection for the cuff connector.

- (i)** If the measurement is performed on the right upper arm, the line should be located on the inside of your elbow. Ensure that your arm is not pressing on the line.


Blood pressure may vary between the right and left arm, which may mean that the measured blood pressure values are different. Always perform the measurement on the same arm.

If the values between the two arms are significantly different, please consult your doctor to determine which arm should be used for the measurement.

Important: The unit may only be operated with the original cuff. The cuff is suitable for an arm circumference of 22 to 36 cm.

A larger cuff for upper-arm circumferences of 35 to 44 cm can be obtained from specialist retailers or from the service address using order number 163.387.

Adopting the correct posture

- Before every measurement, relax for about five minutes. Otherwise deviations can occur.
- You can take the measurement while sitting or lying. Make sure that the cuff is at heart level.
- To take your blood pressure, make sure you are sitting comfortably with your arms and back leaning on something. Do not cross your legs. Place your feet flat on the ground.
- To avoid falsifying the measurement, it is important to remain still during the measurement and not to speak.

Performing the blood pressure measurement

As described above, attach the cuff and adopt the posture in which you want to perform the measurement.

- Press the **START/STOP** button  to start the blood pressure monitor.
All display elements are briefly displayed.


The blood pressure monitor will begin the measurement automatically after 3 seconds. The measurement is taken during the inflation phase.

- You can cancel the measurement at any time by pressing the **START/STOP** button .

As soon as a pulse is found, the pulse symbol  is displayed.

- Systolic pressure, diastolic pressure and pulse readings are displayed. The resting indicator display (see section 6) illuminates in accordance with the positive or negative classification.

- **Er** appears if the measurement could not be performed properly (see section 10 "Error messages/troubleshooting"). Repeat the measurement.

Er

- Now select the desired user memory by pressing the **M1** or **M2** memory buttons. If you do not select a user memory, the measurement is stored in the most recently used user memory. The relevant  or  symbol appears on the display.

- Press the **START/STOP** button ① to switch off the blood pressure monitor. The measurement is then stored in the selected user memory.

If **Bluetooth® data transfer** has been activated, data is transferred after having confirmed the user memory. The blood pressure monitor shows the *Bluetooth®* symbol during data transfer. Press the **START/STOP** button ① once again to switch off the blood pressure monitor.

- ① Please note that you must add the blood pressure monitor in “My devices” in the Beurer “HealthManager” app to enable data transfers. The Beurer “HealthManager” app must be active to allow data transfers. If the current data is not displayed on your smartphone, repeat the transfer as described in Chapter 8.

If you forget to turn off the blood pressure monitor, it will switch off automatically after approximately 3 minutes. In this case too, the value is stored in the selected or most recent user memory and the data is transferred if *Bluetooth®* data transfer has been activated. The *Bluetooth®* symbol appears on the blood pressure monitor display during data transfer.

- Wait at least 5 minutes before taking another measurement!

6. Evaluating results

Cardiac arrhythmia:

This unit can identify potential disruptions of the heart rhythm when measuring and if necessary, indicates this after the measurement with the symbol .

This can be an indicator for arrhythmia. Arrhythmia is a condition in which the heart rhythm is abnormal because of flaws in the bioelectrical system that regulates the heartbeat. The symptoms (skipped or premature heart beats, pulse being slow or too fast) can be caused by factors such as heart disease, age, physical make-up, excess stimulants, stress or lack of sleep. Arrhythmia can only be determined through an examination by your doctor.

If the symbol  is shown on the display after the measurement has been taken, it should be repeated. Please ensure that you rest for 5 minutes beforehand and do not speak or move during the measurement. If the symbol  appears frequently, please consult your doctor.

Self-diagnosis and treatment based on the measurements can be dangerous. Always follow your GP's instructions.

WHO classification:

In accordance with the guidelines/definitions of the World Health Organization and the latest findings, the measurements can be classified and assessed according to the following table.

However, these standard values serve only as a general guideline, as the individual blood pressure varies in different people and different age groups etc.

It is important to consult your doctor regularly for advice. Your doctor will tell you your individual values for normal blood pressure as well as the value above which your blood pressure is classified as dangerous.

The classification on the display and the scale on the unit show which category the recorded blood pressure values fall into. If the values of systole and diastole fall into two different WHO categories (e.g. systole in the 'High normal' category and diastole in the 'Normal' category), the graphical WHO classification on the device always shows the higher category; for the example given this would be 'High normal'.

Blood pressure value category	Systole (in mmHg)	Diastole (in mmHg)	Action
Grade 3: severe hypertension	≥ 180	≥ 110	seek medical attention
Grade 2: moderate hypertension	160–179	100–109	seek medical attention
Grade 1: mild hypertension	140–159	90–99	regular monitoring by doctor
High normal	130–139	85–89	regular monitoring by doctor
Normal	120–129	80–84	self-monitoring
Optimal	< 120	< 80	self-monitoring

Source: WHO, 1999 (World Health Organization)

Resting indicator measurement (using HSD diagnostics)

The most frequent error made when measuring blood pressure is taking the measurement when not at rest (haemodynamic stability), which means that both the systolic and the diastolic blood pressures are incorrect in this case. While measuring the blood pressure, the device automatically determines whether you are at rest or not. If there is no indication that the circulatory system is not sufficiently at rest, the  symbol (haemodynamic stability) lights up green and the measurement can be recorded as a reliable blood pressure at rest value.

GREEN: haemodynamically stable

Measurement of the systolic and diastolic pressure is increased when the circulatory system is sufficiently at rest and is a very reliable indicator of resting blood pressure.

However, if there is an indication that the circulatory system is not sufficiently at rest (haemodynamic instability), the  symbol lights up red.

In this case, the measurement should be repeated after a period of physical and mental rest. The blood pressure measurement must be taken when the patient is physically and mentally rested, as it will be the basis for a diagnosis and regulation of the patient's medical treatment.

RED: lack of haemodynamic stability

It is very probable that the systolic and diastolic blood pressures have not been measured whilst the patient is at rest and the resting blood pressure measurement has therefore been distorted.

Repeat the measurement after a rest and relaxation period of at least five minutes. Go to a sufficiently quiet and comfortable

spot and remain there calmly; close your eyes, breathe deeply and evenly and try to relax.

If the next measurement also shows insufficient stability, you can repeat the measurement after another resting period. If the measurements continue to show some instability, identify these blood pressure measurements as having been taken when the circulatory system had not been sufficiently rested.

In this case, nervousness or inner anxiety may be the cause and this cannot be cured by brief periods of rest. Existing cardiac arrhythmias may also prevent a stable blood pressure measurement.

A lack of resting blood pressure can have various causes, such as physical or mental strain or distraction, speaking or experiencing cardiac arrhythmias during the blood pressure measurement.

In an overwhelming number of cases, the HSD diagnosis will give a very good guide as to whether the circulatory system is rested when taking the measurement. Certain patients suffering from cardiac arrhythmia or chronic mental conditions can remain haemodynamically unstable in the long-term, something which persists even after repeated periods of rest. The accuracy of the resting blood pressure results is reduced in these users. Like any medical measurement method, the precision of the HSD diagnosis is limited and can lead to incorrect results in some cases. The blood pressure measurements taken when the circulatory system was at rest represent particularly reliable results.


7. Displaying and deleting measurements

The results of every successful measurement are stored together with the date and time. The oldest measurement is overwritten in the event of more than 60 measurements.

- To access memory recall mode, the blood pressure monitor must first be started. To do this press the **START/STOP button** ①.
- Within 3 seconds of the full-screen display appearing, select the desired user memory (M1 M2) with the **M1** or **M2** memory button.
 - To view the measurements for user memory M1, press the **M1** memory button.
 - To view the measurements for user memory M2, press the **M2** memory button. Your last measurement will appear on the display.

If *Bluetooth*® is activated (the  symbol appears on the display), measurements are automatically transferred.

Press the **M1/M2** button to cancel data transfer. The device will show average values. The  symbol is no longer shown.


User memory

- ① Press the **M1** memory button if you have selected user memory 1.
Press the **M2** memory button if you have selected user memory 2.

Average values

R flashes on the display.

The average value of all saved measured values in this user memory is displayed.

AM flashes on the display.

The average value of the morning measurements for the last 7 days is displayed (morning: 5.00 a.m. – 9.00 a.m.).

- Press the relevant memory button (**M1** or **M2**).

PM flashes on the display.

The average value of the evening measurements for the last 7 days is displayed (evening: 6.00 p.m. – 8.00 p.m.).

Individual measured values

- When the relevant memory button (**M1** or **M2**) is pressed again, the last individual measurement is displayed (in this example, measurement 03).

- When the relevant memory button (**M1** or **M2**) is pressed again, you can view your individual measurements.
- To switch the device off again, press the START/STOP button **I**.

- i** You can exit the menu at any time by pressing the **START/STOP** button **I**.

Deleting measured values

- To clear the memory of the relevant user memory, you must first select a user memory.
- Start individual measurement access.
- Press and hold the **M1/M2** memory buttons for 5 seconds.

All the values in the current user memory are deleted.

- i** It is not possible to delete individual measurements.

8. Transferring measurements

Transfers using the USB interface

Connect the blood pressure monitor to your PC using the USB cable.

 No data transfer may be launched whilst performing a measurement.

PC is shown on the display.

Begin the data transfer in the “HealthManager” PC software. During the data transfer, an animation is shown on the display. A successful data transfer is displayed as in figure 1. If the data transfer is unsuccessful, an error message appears as in figure 2. In this case, interrupt the PC connection and start the data transfer again.

fig. 1

fig. 2

After 30 seconds of not being in use or if communication with the PC is interrupted, the blood pressure monitor switches itself off automatically.

Bluetooth® Smart transfers

It is also possible to transfer the measured values saved on the device to your smartphone using *Bluetooth®* Smart.

You will need the Beurer „HealthManager“ app for this. The app is available from the App Store.

Proceed as follows to transfer values:

data is transferred if *Bluetooth®* has been activated in the settings menu. The  symbol appears in the top left of the display (see section 4 “Preparing for the measurement”).

Step 1: BM 85

activate *Bluetooth®* on your device (see section 4 “Preparing for the measurement, *Bluetooth®*”).

Step 2: Beurer “HealthManager” app

In the Beurer “HealthManager” app, add the BM 85 under “Settings/My devices”.

Step 3: BM 85

take a measurement.

Step 4: BM 85

data transfer immediately following measurements.

- Select the desired user memory. Start *Bluetooth®* data transfer (chapter 5)

Step 4: BM 85

Data transfer at a later point:

- go to memory mode (chapter 7). Select the desired user memory. *Bluetooth®* data transfer starts automatically.

 The Beurer “HealthManager” app must be active to allow data transfers.

If your smartphone has a protective cover, remove this to ensure that there is no interference during the transfer. Begin the data transfer in the Beurer “HealthManager” app.

9. Cleaning and storing the unit

- Clean the unit and cuff carefully using a slightly damp cloth only.
- Do not use any cleaning agents or solvents.
- Under no circumstances should you hold the unit under water, as this can cause liquid to enter and damage the unit.
- If you store the unit, no heavy objects should be placed on top of it. The cuff line should not be bent sharply.

10. Error messages/trouble-shooting

In case of faults, the E_r _ message appears in the display. Error messages may appear if:

- systolic or diastolic pressure could not be measured ($E_r 1$ or $E_r 2$ appears on the display)
- systolic or diastolic pressure was outside the measurement range (**Hi** or **Lo** appears on the display)
- the cuff is fastened too tightly or loosely ($E_r 3$ or $E_r 4$ appears on the display)
- the pump pressure is higher than 300 mmHg ($E_r 5$ appears on the display)
- pumping up takes longer than 160 seconds ($E_r 6$ appears on the display)
- there is a system or device error ($E_r A$, $E_r U$, $E_r 7$ or $E_r B$ appears on the display)
- the data could not be sent to the PC ($P_C E_r$ appears in the display).

- the data could not be sent via *Bluetooth*® (E_r).

In such cases, repeat the measurement and/or the data transfer. Ensure that the cuff line is correctly attached and that you do not move or speak during the measurement.

Technical alarm – description

Should the recorded blood pressure (systolic or diastolic) lie outside the limits specified in the section “Technical specifications”, the technical alarm will appear on the display indicating either “**Hi**” or “**Lo**”. In such cases, you should seek medical assistance and check the accuracy of your procedure.

The limit values for the technical alarm are factory set and cannot be adjusted or deactivated. These alarm limit values are accorded second priority under the standard IEC 60601-1-8. The technical alarm is a non-locking alarm and must not be reset. The signal shown on the display will disappear automatically after about 8 seconds.

11. Battery

The BM 85 is equipped with a lithium ion battery (3.7V/400 mAh). If the  symbol flashes, you must charge the battery for a minimum of two hours using the enclosed cable. Recharge the battery to between 50 and 75 % of its capacity at least twice a year to achieve a maximum battery service life.

Device charged (75% - 100%)

Device charged (50% - 75%)

Device charged (25% - 50%)

	< 25%
	You can take a maximum of 10 measurements (flashes)

The *Bluetooth*® word mark and logos are registered trademarks owned by *Bluetooth*® SIG, Inc. and any use of such marks by Beurer GmbH is under license. Other trademarks and trade names are those of their respective owners.

12. Technical specifications

Model no.	BM 85
Measurement method	Oscillometric, non-invasive blood pressure measurement on the upper arm
Measurement range	Cuff pressure 0–300 mmHg, systolic 60–260 mmHg, diastolic 40–199 mmHg, pulse 40–180 beats/minute
Display accuracy	Systolic ± 3 mmHg, diastolic ± 3 mmHg, pulse ± 5 % of the value shown
Measurement inaccuracy	Max. permissible standard deviation according to clinical testing: systolic 8 mmHg/ diastolic 8 mmHg
Memory	2 x 60 memory spaces
Dimensions	L 180 mm x W 99 mm x H 40 mm
Weight	Approximately 317 g (without cuff)

Cuff size	22 to 36 cm
Permissible operating conditions	+10 °C to +40 °C, ≤ 90 % relative humidity (non-condensing)
Permissible storage conditions	-20 °C to +55 °C, ≤ 90 % relative humidity, 800–1050 hPa ambient pressure
Power supply	DC 5V --- 600mA lithium ion battery 3.7 V/400mAh
Battery life	For approx. 50 measurements, depending on levels of blood pressure and pump pressure
Accessories	Cuff, cuff holder, instructions for use, USB mains part, USB cable, storage pouch
Classification	Internal supply, IPX0, no AP or APG, continuous operation, application part type BF
Data transfer via <i>Bluetooth</i> ® wireless technology	The blood pressure monitor uses <i>Bluetooth</i> ® Smart (low energy), 2.4 GHz frequency band, Compatible with Bluetooth 4.0 smartphones/tablets
	List of supported smartphones/tablets

Technical information is subject to change without notification to allow for updates.

- This device complies with European Standard EN60601-1-2 and is subject to particular precautions with regard to electromagnetic compatibility. Please note that portable and mobile HF communication systems may interfere with this unit. More details can be requested from the stated Customer Services address or found at the end of the instructions for use.
- This device corresponds to the EU Medical Devices Directive 93/42/EC, the German Medical Devices Act (Medizinproduktegesetz) and the standards EN 1060-1 (non-invasive sphygmomanometers, Part 1: General requirements), EN 1060-3 (non-invasive sphygmomanometers, Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems) and IEC 80601-2-30 (Medical electrical equipment – Part 2–30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers).
- The accuracy of this blood pressure monitor has been carefully checked and developed with regard to a long useful life. If using the device for commercial medical purposes, it must be regularly tested for accuracy by appropriate means. Precise instructions for checking accuracy may be requested from the service address.
- We hereby guarantee that this product complies with the European R&TTE Directive 1999/5/EC. Please contact the specified service address to obtain further information, such as the CE Declaration of Conformity.

13. Mains part

Model No.	OH-1048A0500600U2 VDE
Input	100–240V AC, 50–60 Hz; 120-60 mA
Output	5V DC, 600 mA, in conjunction with Beurer blood pressure monitors only
Manufacturer	Dongguan Oriental Hero Ele. Co., LTD.
Protection	The device is double protected and has a primary-side cutout switch which disconnects the device from the mains in case of malfunction.
	Polarity of the the DC voltage connection
	Insulated / protection class 2
Housing and protective covers	The housing of the mains part protects users from touching live parts or parts that could be live (for example with their fingers, or with a needle or checking hook). The user must not touch the patient and the output connector of the AC mains part at the same time.

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Table 1

For all ME EQUIPMENT and ME SYSTEMS

Guidance and manufacture's declaration – electromagnetic emission

The BM 85 is intended for use in the electromagnetic environment specified below. The customer or the user of the BM 85 should assure that it is used in such an environment.

Emission test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The BM 85 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The BM 85 is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 2**For all ME EQUIPMENT and ME SYSTEMS***Guidance and manufacture's declaration – electromagnetic immunity*

The BM 85 is intended for use in the electromagnetic environment specified below. The customer or the user of the BM 85 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	<5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles <5 % U_T (>95 % dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the BM 85 requires continued operation during power mains interruptions, it is recommended that the BM 85 be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			

Table 3
For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Guidance and manufacture's declaration – electromagnetic immunity

The BM 85 is intended for use in the electromagnetic environment specified below. The customer or the user of the BM 85 should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	<p>Portable and mobile RF communications equipment should be used no closer to any part of the BM 85, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	<p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,^a should be less than the compliance level in each frequency range.^b</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol: </p>

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BM 85 is used exceeds the applicable RF compliance level above, the BM 85 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BM 85.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.

Table 4
For ME EQUIPMENT and ME SYSTEMS that are not LIFE-SUPPORTING

Recommended separation distances between portable and mobile RF communications equipment and the BM 85

The BM 85 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BM 85 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BM 85 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.